

# PureBeam™ delivery system

The PureBeam delivery system is the preferred tool of minimally invasive surgery specialists to treat soft tissue disorders in the field of Neuro-, Spine-, and ENT surgery.

**SureGuide™**  
Laser energy delivery system propagates laser light in TEM<sub>00</sub> mode through a thin flexible and robust titanium alloy jacketed cable.

**SureKit™**  
Laser turret assembly that enables the precise launch of a TEM<sub>00</sub> mode beam from various CO<sub>2</sub> laser systems into the SureGuide.

**SureLase™**  
CO<sub>2</sub> laser system for delivery of stable TEM<sub>00</sub> mode laser power up to 30 Watts in continuous, gated or super pulse mode.

**SureBlade™**  
Ultra-sharp, durable and sterilizable diamond laser scalpel for precision surgeries with concurrent coagulation.

**SureScan™**  
High velocity laser scanner system for precise computer controlled ablation of skin, scars and adhesions.

**Added security —**  
the cumbersome articulated arm is replaced with a waveguide cable that will not go out of alignment.



 **CLINICON™**  
PRECISION FOCUSED™

[purebeam@clinicon.com](mailto:purebeam@clinicon.com)  
[www.clinicon.com](http://www.clinicon.com)

# PureBeam™ delivery system

The Clinicon SureTouch universal handpiece and SureProbe tips are used to deliver CO<sub>2</sub> laser energy for the incision, excision, vaporization, ablation, coagulation or cauterization of soft tissue. Examples of some procedures in ENT, Neuro, and Spine are:

## ENT

Pulmonary Lesions  
Vocal Cord Dysplasia & Carcinoma  
Otosclerosis (Stapedectomy)  
Turbinectomy  
Bronchoectomy  
Cordectomy  
Polyps  
Larynx Leukoplakia  
Arytenoidectomy

Papillomata  
Phenoidal Pituitary Tumor Removal  
Lymphangioma  
Hemiglossectomy  
Hypophysectomy  
Tuberoseptum Synechia  
Acoustic Neuroma  
Laser Assisted Myringotomy (LAM)  
Laser Assisted Uvulo-Palatoplasty (LAUP)


In-excision Biopsy  
Removal of Benign & Malignant Lesions  
Frenectomy  
Gingivectomy  
Gingivoplasty  
Tuberosity Reduction  
Operculectomy  
Aphthous Ulcer Therapy  
Coagulation of Graft Donor Sites

## NEURO

Glioma  
Glioblastoma  
Medulloblastoma  
Astrocytoma  
Metastatic Lesions  
Craniopharyngioma  
Meningioma  
Acoustic Neuroma  
Pituitary Tumors

## SPINAL

Benign and Malignant Tumors  
Astrocytomas  
Lymphomas  
Meningiomas  
Chordomas  
Sarcomas



## PART# DESCRIPTION

- 3590** SureGuide™ 275cm Fiber Cable, Mechanical, Packaged Non-Sterile Assembly
- 3156** SureGuide™ 275cm Fiber Cable, Electrical, Packaged Non-Sterile Assembly
- 4033** SureGuide™ 200cm Fiber Cable, Mechanical, Packaged Non-Sterile Assembly
- 4034** SureGuide™ 200cm Fiber Cable, Electrical, Packaged Non-Sterile Assembly
- 3730** SureTouch™ 1.5mm Universal Handpiece, Packaged Non-Sterile Assembly



- 3901** SureTip™ FREE-BEAM, Straight Insert, Packaged Sterile in Box of 10

- 4050** SureKit™ ARC, Laser Turret Adapter, Packaged Assembly

- 2681** SureKit™ LUXAR, Laser Turret Adapter, Packaged Assembly

- 4052** SureKit™ SHARPLAN, Laser Turret Adapter, Packaged Assembly

- 4060** SureLase™ 30 Watt CO<sub>2</sub> Laser System




- 3566** SureScan™ Laser Scanner, Packaged Assembly



- 2572** SureBlade™ LANCE, 6mm Diamond Laser Scalpel, Packaged Non-Sterile Assembly

- 2986** SureBlade™ BI-FACET, 8mm Diamond Laser Scalpel, Packaged Non-Sterile Assembly



SureProbe™ is available in a dispensing box that holds ten individually-sterilized probe pouches

FOR MORE INFORMATION OR TO PLACE AN ORDER CONTACT  
1-800-CLINICON OR VISIT OUR WEBSITE AT WWW.CLINICON.COM



**Clinicon Corporation USA**  
3025-A Industry Street  
Oceanside, CA 92054  
T 800-CLINICON  
T 760-439-1700  
F 760-439-1798  
purebeam@clinicon.com  
www.clinicon.com