

PureBeam™ delivery system

The PureBeam delivery system is the preferred tool of minimally invasive surgery specialists to treat soft tissue disorders in the field of Neuro-, Spine-, and ENT surgery.

SureGuide™
Laser energy delivery system propagates laser light in TEM₀₀ mode through a thin flexible and robust titanium alloy jacketed cable.

SureKit™
Laser turret assembly that enables the precise launch of a TEM₀₀ mode beam from various CO₂ laser systems into the SureGuide.

SureLase™
CO₂ laser system for delivery of stable TEM₀₀ mode laser power up to 30 Watts in continuous, gated or super pulse mode.

SureBlade™
Ultra-sharp, durable and sterilizable diamond laser scalpel for precision surgeries with concurrent coagulation.

SureScan™
High velocity laser scanner system for precise computer controlled ablation of skin, scars and adhesions.

Added security —
the cumbersome articulated arm is replaced with a waveguide cable that will not go out of alignment.



 **CLINICON™**
PRECISION FOCUSED™

purebeam@clinicon.com
www.clinicon.com

PureBeam™ delivery system

The Clinicon SureTouch universal handpiece and SureProbe tips are used to deliver CO₂ laser energy for the incision, excision, vaporization, ablation, coagulation or cauterization of soft tissue. Examples of some procedures in ENT, Neuro, and Spine are:

ENT

Pulmonary Lesions
Vocal Cord Dysplasia & Carcinoma
Otosclerosis (Stapedectomy)
Turbinectomy
Bronchoectomy
Cordectomy
Polyps
Larynx Leukoplakia
Arytenoidectomy

Papillomata
Pheoidal Pituitary Tumor Removal
Lymphangioma
Hemiglossectomy
Hypophysectomy
Tuberoseptum Synechia
Acoustic Neuroma
Laser Assisted Myringotomy (LAM)
Laser Assisted Uvulo-Palatoplasty (LAUP)


In-excision Biopsy
Removal of Benign & Malignant Lesions
Frenectomy
Gingivectomy
Gingivoplasty
Tuberosity Reduction
Operculectomy
Aphthous Ulcer Therapy
Coagulation of Graft Donor Sites

NEURO

Glioma
Glioblastoma
Medulloblastoma
Astrocytoma
Metastatic Lesions
Craniopharyngioma
Meningioma
Acoustic Neuroma
Pituitary Tumors

SPINAL

Benign and Malignant Tumors
Astrocytomas
Lymphomas
Meningiomas
Chordomas
Sarcomas



PART# DESCRIPTION

- 3590** SureGuide™ 275cm Fiber Cable, Mechanical, Packaged Non-Sterile Assembly
- 3156** SureGuide™ 275cm Fiber Cable, Electrical, Packaged Non-Sterile Assembly
- 4033** SureGuide™ 200cm Fiber Cable, Mechanical, Packaged Non-Sterile Assembly
- 4034** SureGuide™ 200cm Fiber Cable, Electrical, Packaged Non-Sterile Assembly
- 3730** SureTouch™ 1.5mm Universal Handpiece, Packaged Non-Sterile Assembly



- 3901** SureTip™ FREE-BEAM, Straight Insert, Packaged Sterile in Box of 10

- 4050** SureKit™ ARC, Laser Turret Adapter, Packaged Assembly

- 2681** SureKit™ LUXAR, Laser Turret Adapter, Packaged Assembly

- 4052** SureKit™ SHARPLAN, Laser Turret Adapter, Packaged Assembly

- 4060** SureLase™ 30 Watt CO₂ Laser System




- 3566** SureScan™ Laser Scanner, Packaged Assembly



- 2572** SureBlade™ LANCE, 6mm Diamond Laser Scalpel, Packaged Non-Sterile Assembly

- 2986** SureBlade™ BI-FACET, 8mm Diamond Laser Scalpel, Packaged Non-Sterile Assembly



SureProbe™ is available in a dispensing box that holds ten individually-sterilized probe pouches

FOR MORE INFORMATION OR TO PLACE AN ORDER CONTACT
1-800-CLINICON OR VISIT OUR WEBSITE AT WWW.CLINICON.COM



Clinicon Corporation USA
3025-A Industry Street
Oceanside, CA 92054
T 800-CLINICON
T 760-439-1700
F 760-439-1798
purebeam@clinicon.com
www.clinicon.com